

User data report

Daily report from Jan 26, 6am – Jan 27, 6ams

Overview Keexxx & VivoSmart 5

Keexxx



Collected data: 13 hours

- Night: 8 hours
- Day: 6 hours
- Rest: 1 hours
- Activity: 1 hours


Garmin watch



Collected data: 13 hours

- Night: 8 hours
- Day: 8 hours
- Rest: 2 hours
- Activity: 3 hours

Duration of data collection



Health overview

Cardiology

- Tachycardia
- Bradycardia
- QTc Length anomaly
- Heart Rate (HR) XXms
- HR max 119bpm
- HR min 51bpm
- Sdnn (HRV) XXms

Respiratory

- Tachypnea
- Bradypnea
- Breath Rate XXms
- BRV XXms

Spo2

- Average Spo2
- 91%**
- 84%** lowest
- <70%
 - 70-79%
 - 80-89%
 - 90-100%

Sleep

- Sleep score
- 80/100**
- Deep: 10%
 - Light: 40%
 - Paradoxal: 10%
 - Awake: 40%

Stress

- Stress score
- 37**
- Rest 2h20 (12%)
 - Low 1h53 (10%)
 - Medium 3h48 (20%)
 - High 10h05 (55%)

Steps

- Number of steps
- 8230**
- 450** steps on average per day
- 3,9km** distance

Body Battery

- High
- 100/100**
- Low
- 11/100**

Calories


- Number of calories
- 7,100 kcals**
- + mean number of calories per day if weekly or monthly

User data report


Daily report from Jan 26, 6am – Jan 27, 6ams

Keexxx Cardiology data


Beats per minute



Time interval between beats



QT interval length - Corrected with Framingham formula



User data report

Daily report from Jan 26, 6am – Jan 27, 6ams

Period & Activity Definition

-  **Night:** From 0am to 6am
-  **Day:** From 6am to 0pm
-  **Rest:** Resting or low activity during the day
-  **Activity:** Moderate or high motion detection

Cardiac Indicators Definitions

- BPM**
Number of Heart Beats Per Minute
- HR**
Time in milliseconds between successive QRS complexes (RR intervals)
- Sdnn**
Standard deviation of RR intervals on 5min segments
- QTc interval**
Time between Q and T waves in milliseconds, normalised by Framingham formula:
 $1000 * (QT/1000 + 0.154 * (1-RR))$
- pNN50**
pNN50 value quantifies the percentage of adjacent NN intervals that differ by more then 50ms

Respiratory Indicators Definitions

- RPM**
Number of Respiratory Cycles Per Minute
- BR**
Breath Rate. Time in seconds between successive respiratory cycles
- BRV**
Breath Rate Variability. Standard deviation of breath rate

Alert Notifications

-  A red check box is used to notify if there are any alerts

Data Quality

Data Quality Index Percentage of exploitable data after data cleaning, averaged over all days of measure